

## IS YOUR HEARING AID IN, ON, AND WORKING?

EVERY DAY		
<b>TURN YOUR HEARING AID ON EACH MORNING</b>		<b>DAY:</b> Turn <b>ON</b> by closing the battery door
<b>WEAR YOUR HEARING AID EVERY DAY – GET INTO A ROUTINE</b>		Put your ear <b>mould in correctly</b> . Is it comfortable? It should not whistle or buzz.
EVERY NIGHT		
<b>SWITCH OFF YOUR HEARING AID EACH NIGHT</b>		<b>NIGHT:</b> Turn <b>OFF</b> by opening the battery door
<b>CLEAN YOUR EAR MOULD EVERY NIGHT</b>	or 	Clean your ear mould with an unscented <b>baby wipe</b> or <b>clean damp cloth</b> to stop it <b>getting blocked</b> .
EVERY WEEK		
<b>CHANGE YOUR BATTERIES WEEKLY (they last 7-10 days)</b>		Remove the label from the new battery and <b>put the battery in + side up</b>
EVERY 6 MONTHS		
<b>REPLACE THE TUBING</b>		<b>Contact Audiology</b> to get the tubing changed.

Shropshire NHS Audiology Telephone: 01743 261482 [sath.audiology@nhs.net](mailto:sath.audiology@nhs.net)

For Information: [www.sath.nhs.uk/services/audiology](http://www.sath.nhs.uk/services/audiology) and [www.C2hearOnline.com](http://www.C2hearOnline.com)

## HEARING AID DETAILS FOR: \_\_\_\_\_

I have a hearing loss and wearing hearing aid/s Date: \_\_\_\_\_

### Hearing Aid Batteries

<b>I can get my hearing aid batteries from:</b>	My local GP practice:  Other place e.g. Hospital:  By ringing Audiology on 01743 261482/Send brown book in the post with a request for batteries  Email: <a href="mailto:sath.audiology@nhs.net">sath.audiology@nhs.net</a>  Self-request via electronic request on Audiology website <a href="http://www.sath.nhs.uk/audiology">www.sath.nhs.uk/audiology</a>
<b>What Type of Batteries?</b>  <b>Where do I keep my batteries?</b>	Battery Code: P10 P13 P312 P672 (circle) Battery Label Colour: Yellow Orange Brown Blue

### Hearing Aid Care and Maintenance

<b>My local NHS Audiology clinic is at:</b>	
<b>My hearing aid details:</b>	Left: _____ Right: _____
<b>Tubing Type or code (Dome if appropriate):</b>	
<b>Is the volume control active?</b>	Yes _____ No _____ (Circle) Details if yes: _____
<b>Hearing aid programme details:</b>	
<b>Where are the hearing aids kept when not in use?</b>	Colour of hearing aid box: _____ Where kept: _____