

Family/Friends/Carers Supporting Hearing Aid Use

Hear better, strain less, stimulate the brain with sound, keep socially active

| EVERY DAY | | |
|--|-------------|--|
| DAYTIME : Turn ON by closing the battery door | | Check the hearing aid is working by switching it on, cupping your hands around the mould and aid. Listen for a whistle/buzz. |
| PROMOTE WEARING HEARING AIDS DAILY AS ROUTING TO ENSURE THE BRAIN IS STIMULATED BY EVERYDAY SOUNDS | Hearing ald | Help to put the ear mould in correctly – you can gently pull the ear lobe down. Wearing the hearing aid/s should be part of a daily living routine eg after washing in the morning. |
| EVERY NIGHT | | |
| NIGHT : Turn OFF by opening the battery door | OFF Click | Keep the hearing aid box in the same place at night in a visible place eg bedside table. It could have left/right hearing aid written on the box. Red dot for Right ear. |
| CLEAN THE EAR MOULD EVERY NIGHT | OR So | Clean the ear mould with a baby wipe or clean damp cloth to stop it getting blocked with wax or dry skin and keep it working well. |
| EVERY WEEK | | |
| CHANGE BATTERIES WEEKLY eg every Sunday (They can last 7 – 10 days and beep when running out) | | Remove the label from the new battery and put the battery in + side up. Check each month how many batteries are left to make sure they don't run out and ring/email audiology for more or get more from your GP surgery.Online ordering too. |
| EVERY 6 MONTHS | | |
| REPLACE THE TUBING | | Put reminders on the calendar every 6 months. Contact Audiology to get the tubing changed |
| Shropshire NHS Audiology Telephone: 01743 261482 <u>sath.audiology@nhs.net</u> For Information: <u>www.sath.nhs.uk/services/audiology</u> | | |
| www.C2Hearonline.com: you Tube Clins on bearing aid care and maintenance | | |

www.C2Hearonline.com: you Tube Clips on hearing aid care and maintenance

Patient Advice and Liaison Services: 01952 282888/01743261691. If you require this information in large print or another language please contact audiology. HA Poster Carer instructions v8