

IS YOUR HEARING AID IN, ON, AND WORKING?

EVERY DAY		
TURN YOUR HEARING AID ON EACH MORNING		DAY: Turn ON by closing the battery door
WEAR YOUR HEARING AID EVERY DAY – GET INTO A ROUTINE		Put your ear mould in correctly . Is it comfortable? It should not whistle or buzz.
EVERY NIGHT		
SWITCH OFF YOUR HEARING AID EACH NIGHT		NIGHT: Turn OFF by opening the battery door
CLEAN YOUR EAR MOULD EVERY NIGHT	or 	Clean your ear mould with an unscented baby wipe or clean damp cloth to stop it getting blocked .
EVERY WEEK		
CHANGE YOUR BATTERIES WEEKLY (they last 7-10 days)		Remove the label from the new battery and put the battery in + side up
At LEAST EVERY 6 MONTHS		
REPLACE THE TUBING		Contact Audiology to get the tubing changed.

Shropshire NHS Audiology Telephone: 01743 261482 sath.audiology@nhs.net
For Information: www.sath.nhs.uk/services/audiology & www.C2hearOnline.com